

Automated Closed Cup Flash Point Analyzers

HFP 36X Series



HERZOG
by **PAC**

Proven and Reliable Flash Point Testing

- Pensky Martens, TAG and Abel tests with a single instrument
 - Easy and safe operation
 - Extended program library and results database
- State-of-the-art calibration and diagnostics features

PROVEN, RELIABLE AND COST-EFFECTIVE FLASH POINT TESTING

The series of HFP Flash Point Analyzers is well known by their proven reliability, superior design and quality. The HFP series has been designed to save you time, money and bench space. The HFP automatically determines flash point in strict compliance with the appropriate test method. In addition to standard test methods, the HFP may be programmed with as many as 19 user-defined test protocols and to test for samples with an unknown flash point.

To increase test productivity on products having elevated flash point, there is a capability to start with fast heating rate. The instrument automatically switches to the standard heating rate at appropriate temperature to assure accurate result. When flash point is detected, results are automatically corrected for standard barometric pressure and displayed on a bright, easy to read HFP's graphics readout and can also be output automatically on an optional printer or computer system. Herzog token ring system allows easy resource sharing and facilitate LIMS connection.

The series of interchangeable test modules allow quick, easy switches from one test method to another. The HFP automatically recognizes the type of module installed, and you can begin another standard method testing in minutes.



HFP Series Automatic Flash Point Analyzers are designed to simplify the testing of hydrocarbon samples using the Pensky-Martens or other closed cup test methods.

APPLICATION RANGE

Volatility

- Petroleum products
- Biodiesels
- Solvents
- Chemicals
- Fluxed bitumen
- Food and Beverage

STANDARD TEST METHODS

Pensky Martens:

- ASTM D93 A, B, C
- EN ISO 2719 A, B, C
- DIN 51758
- IP 34 A, B, C
- JIS K 2265

Tag:

- ASTM D56

Abel:

- IP 170
- ISO 13736

RELIABLE FLASH POINT ANALYSES



- Strictly respects all requirements of test methods
- Proven thermal detection precisely captures flash point eliminating interference from water or silicone
- Tests for samples with unknown flash points
- Pt-100 temperature sensor duplicates response time of a mercury thermometer. Results are automatically corrected for standard barometric pressure and display on a bright, easy-to-read graphic screen
- Built-in air cooler quickly cools heating block in preparation for next test; optional cooling block enables sub-ambient Pensky-Martens testing
- Special program assures precise testing of viscous samples, like fluxed bitumens
- Integral database stores up to 300 test results with statistical analysis functions

SAFE, DEPENDABLE OPERATION

- Overheating protection system with automatic heating shut-off ensures safety during test; automatic fire detection system provides audible alarm with external alarm connection available
- Relights flame, if necessary, during operation and suppresses gas source at end of test
- Extensive diagnostics allow user to confirm and manually operate all analyzer functions; electronic temperature measurement circuit continually self-calibrates
- Quality construction and reliable operation backed by a limited parts and service warranty
- Expert sales & service provided by PAC's worldwide network of factory-trained, authorized representatives

OPERATION FLEXIBILITY

- **Ignition:** HFP Analyzers are available with either an electric or gas ignition system. Units equipped with a special circuit to reduce and monitor electric ignitor wear
- **Flash detection:** The HFP is equipped with a thermal detection system, eliminating interference from water or silicone-containing samples
- **Fire safety:** A built-in fire sensor alerts you instantly to flames outside the flash cup. A potential-free alarm relay contact is also available to link the HFP to a fire suppression or remote alarm system
- **Cooling:** Very efficient built-in air fan in Pensky-Martens module enables fast cooling the heating block to a pre-programmed temperature once test is completed. Required external cooling, TAG and Abel test modules are equipped with fast coupling connectors that minimize installation time to a few minutes
- **User interface:** Longlife full color touch screen display with screen saver and replaceable protection film. Easy navigation through operational procedures owing to modern 'look and feel' machine interface. Pop-up online help is always at your hand to inquire
- **Versatility:** In addition to performing standard Pensky-Martens, TAG and Abel testing, this versatile analyzer is easily programmed with customized test protocols



The HFP's touch screen menu lets you quickly initiate standard flash point testing, recall test methods, change test parameters, run the system diagnostics, read messages, and go through calibrations.

Flash point determinations may be run using the standard test method or customized test protocols. Standard methods include Pensky-Martens (ASTM D93 A,B,C), TAG Closed Cup (ASTM D56) and Abel (IP 170).

Easy to view the last 20 tested samples—a helpful tool if the same samples are run on a daily basis. Select the sample name, which automatically recalls method and expected flash point, then press Start. For unknown samples, mark Pretest to prevent dangerous situation with wrong expected flash point entry.

The HFP automatically corrects the measured flash point temperature for standard barometric pressure. Also, test status, sample and block temperature are indicated.

The HFP incorporates comprehensive calibration features. Calibration is time tagged with lock control. Sample probe multi-point calibration table allows separate electronic and probe offsets tracking for up to 12 temperature points in the range of -50 to 400°C. Calibration History log file can be viewed and printed.



U.S.A.

PAC, LP | 8824 Fallbrook Drive | Houston, Texas 77064
T: +1 800.444.TEST | O: +1 281.940.1803 | F: +1 281.580.0719
sales.usa@pacpl.com | service.usa@pacpl.com

FRANCE

BP 70285 | Verson | 14653 CARPIQUET Cedex
T: +33 (0) 231 264 300 | F: +33 (0) 321 266 293
sales.france@pacpl.com | service.france@pacpl.com

GERMANY

Badstrasse 3-5 | P.O.Box 1241 | D-97912 Lauda-Königshofen,
T: +49 9343 6400 | F: +49 9343 640 101
sales.germany@pacpl.com | service.germany@pacpl.com

SINGAPORE

61 Science Park Road | #03-09/10 The Galen
Singapore Science Park III | Singapore 117525
T: +65 6412 0890 | F: +65 6412 0899
sales.singapore@pacpl.com | service.singapore@pacpl.com

The NETHERLANDS

P.O.Box 10.054 | 3004 AB Rotterdam
Innsbruckweg 35 | 3047 AG Rotterdam
T: +31 10 462 4811 | F: +31 10 462 6330
sales.netherlands@pacpl.com | service.netherlands@pacpl.com

RUSSIA

Shabolovka Street | 34, Bldg. 2 | 115419 Moscow
T: +7 495 617 10 86 | F: +7 495 913 97 65
sales.russia@pacpl.com | service.russia@pacpl.com

CHINA

Room 1003, Sunjoy Mansion | No. 6 RiTan Rd.
Chao Yang District | Beijing 100020
T: +86 10 650 72236 | F: +86 10 650 72454
sales.china@pacpl.com | service.china@pacpl.com

INDIA

1508 | Dev Corpora | Pokhran Road No.1
Eastern Express Highway | Thane (W) - 400 601
T: +91-22-6700 4848 | F: +91-22-4228 4950
sales.india@pacpl.com | service.india@pacpl.com

MIDDLE EAST

A1 Quds Street, A1 Tawar road | LIU#H13 Dubai Airport Freezone
Near Dubai Airport (terminal 2) | P.O.Box #54781 | Dubai, UAE
T: +971 04 2947 995 | F: +971 04 2395 465
sales.middleeast@pacpl.com | service.middleeast@pacpl.com

SOUTH KOREA

#621 World Vision Building | 24-2, Youido-dong
Seoul 150-010
T: +82 2785 3900 | F: +82 2785 3977
sales.southkorea@pacpl.com | service.southkorea@pacpl.com

THAILAND

26th Floor, M. Thai Tower | All Seasons Place
87 Wireless Road | Lumpini, Phatumwan | Bangkok 10330
T: +66 2627 9410 | F: +662627 9401
sales.thailand@pacpl.com | service.thailand@pacpl.com

PAC Authorized Representatives are also located in most countries worldwide. For more information visit www.pacpl.com

HERZOG BY PAC

Herzog, originally established in 1937, is a long-established comprehensive line of laboratory instruments for testing distillation, flash point, vapor pressure, bitumen testing, cold flow properties, viscosity and other physical properties of materials.

SPECIFICATIONS

General Information	
Ordering Information	HFP 36x Series Base Unit includes touch-screen control display, electric ignition (standard) or gas ignition (optional), and all connections necessary to support HFP 36x Series interchangeable test modules. Base Unit ships with one test module. HFP 360 Pensky-Martens Flash Point Analyzer : ASTM D93 A,B,C; ISO 2719 A,B,C; DIN 22719 A,B,C, DIN 51758; IP 34 A,B, C; JIS K 2265 HFP 362 TAG Flash Point Analyzer: ASTM D56 HFP 364 Abel Flash Point Analyzer: IP 170; ISO 13736
Operation	
Detection	Thermal; eliminates interference from water or silicone-containing samples
Temperature Range	Per appropriate method or user-defined
Heating Range	Per appropriate method or user-defined; automatic proven algorithm controls heat ignition frequency per block and sample temperatures
Sample Stirring	Per appropriate method or user-defined
Ignition Frequency	Per appropriate method or user-defined
Barometric Pressure	Built-in barometric pressure gauge; flash temperature data automatically corrected for barometric pressure
Interface	Large 6", full graphic, 128 color TFT touch screen; operation of instrument is controlled via touch-screen; connection of external keyboard also possible
Cooling	Built-in ventilation; external cooling block (Pensky-Martens only) or cooling connection (TAG and Abel) available for sub-ambient testing, see Options and Accessories below
Password Security	Multi-level password protection
Documentation	
	Local on-screen display (on Base Unit); Base Unit locally stores up to 300 test results and up to 20 standard pre-programmed and user-defined methods; parallel port for optional printer; RS 232 serial port
Diagnostics & Calibration	
	Built-in comprehensive diagnostics program checks every key component and assembly; routine calibration procedures, including self-calibration of sample temperature measuring system
Utility Requirements	
	115 or 230 VAC, 50/60 Hz, 1100 watts
Dimensions & Weight	
	25 cm W x 56 cm D x 51 cm H (24 kg) 9.85" W x 22" D x 20.1" H (53 lbs)
Options and Accessories	
Small Volume (20ml) Testing (HFP 360 Pensky Martens only)	Sample Cup for small volume (20ml) testing P/N 302-029 Sample Cup Cover, electric ignition only P/N 103-066
Sub-Ambient Testing (Flash Point below 40°C)	HFP 360 Pensky Martens: cooling block test cup -cooled with dry ice P/N 302-027 -cooled with built-in coil, for use with external chiller P/N 302-028 HFP 362 TAG & HFP 364 Abel: modules are available with quick-coupling cooling connections. For P/N see standard quotation
External Keyboard	
Certified Reference Material	Wide range of CRM flash point performance material is available from PAC

Continuing research and development may result in specifications or appearance changes at any time